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EN

Indication for use:

Preventing gait impairments (foot drop) in case of the peroneal nerve neuropathy, as well as peroneal nerve paralysis of various aetiology.

Contraindications:

Circulatory disorders and sensory disturbances within the area of using the orthosis (with or without skin lesions). Lymphoedema.

NOTE:

The orthosis should always be used according to the recommendations stated by a physician or physiotherapist.

As far as the first use of orthosis is concerned, it should be put on and adjusted by a qualified health care specialist. Do not put the orthosis directly on bare skin. The orthosis should be put on a cotton sock. Before using the orthosis, it is crucial to carefully read the instructions for use.

If you experience any discomfort when using the orthosis, there is a swelling or skin irritation, please immediately contact a physician or physiotherapist.

FOOT-RISE is a light, dynamic drop foot orthosis and an orthopaedic appliance supporting the dorsiflexion of the foot when this is required to maintain the proper physiology of gait. The basic version comprises a strap encompassing the ankle joint, a transparent insert and a band enabling to use the orthosis without shoes. The strap is made of a space fabric and a flexible foam, which prevents skin irritations. The insert is placed between the shoe tongue and the shoelaces. The position of the insert determines the degree of support as far as dorsiflexion of the foot is concerned. These elements are joined together by means of a strong, flexible strap with a user-friendly fastening system. This particular construction enables lifting the forefoot (supporting the dorsiflexion of the foot) during the phase of moving the limb, which allows to achieve a gait close to the physiological one.

Putting the orthosis on:

The orthosis may be used with shoes (insert), as well as directly on the foot (band encompassing the forefoot).

Detach the insert from the ankle joint strap. Place the ankle joint strap around ankles, bearing in mind that the narrower strap has to be placed at the bottom and the clip for the insert has to be at the front. Fasten the strap tightly, but loosely enough to ensure normal blood flow and to avoid abrasions in the ankle area. Place the insert between the shoe tongue and shoelaces and pull the flexible band between the two rows of shoelaces. The lower you place the insert

(closer to the front of the shoe) the stronger its action (greater dorsiflexion of the foot).

Place the ankle joint strap as in the version anticipated for shoes - around your ankles, with the narrower strap placed at the bottom, and the clip for the insert placed at the front. Place the soft frame for the forefoot (band) immediately behind the mid-foot, and make sure that the clip is placed at the top. Connect the band with the ankle joint strap with the use of the clip.

Size:

S	18 - 21 cm
M	22 - 26 cm
L	27 - 33 cm

Maintenance:

Orthosis can be hand washed (30°) with the use of a gentle laundry detergents. Dry away from artificial sources of heat. Please check the orthosis carefully before using it after washing. Failure to follow proper personal hygiene, using inappropriate laundry detergents may result in abrasions, as well as other conditions related with excessive sweating and the development of bacterial flora.

NOTE:

THE FINAL DECISION CONCERNING THE USE AND SELECTION OF THE DEVICE IS MADE BY A SPECIALISED PHYSICIAN.