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ENG

ASSINGMENT:

Restor is a long wrap-around brace of multifunctional, therapeutic use. It effectively supports treatment and rehabilitation for injuries and dysfunctions around the knee joint, in cases not requiring immobilisation.

INDICATIONS:

Conditions which require immobilization or control of knee range of movement after surgeries or injuries. Treatment and further treatment of knee joint and soft tissue injuries.

**CONTRAINDICATIONS:**

- Conditions which require cast immobilization.
- Long-healing wounds and ulcers.
- Significant bone deformations (brace cannot be adjusted to the joint).

PRODUCT ADJUSTMENT:

1. Wrap and fasten the foam wraps around the leg using the side bar as a gauge for positioning.
2. Position and press into place the side bars. Reposition the wraps if necessary to align with straps. The hinges should be at the mid-patella and towards the back of the knee on sides. For illustration purposes the straps are not shown in Step 2. Do not remove the straps for application.
3. Secure and fasten the straps closest to the knee first, then the remaining straps.
4. A properly applied post-op knee brace. The hinges allows controlled range of motion with easy adjust pin stops.

5. Extension at: -15°, 0°, 15°, 30°, 45°, 60°, 75°, 90° Flexion at: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°. Provides post-operative immobilization and range of motion needed after ACL, PCL, MCL and LCL surgeries.

**ATTENTION!**

The application instructions listed above are intended for the use of qualified professionals only. It is not recommended for patients to attempt to apply or adjust this device unless otherwise instructed to do by an orthopedic professional. Any questions regarding further adjustment or frequency/ duration of wear should be directed to your orthopedic professional. as they know the individual condition for which this device was prescribed.

SIZE: UNIVERSAL

CLEANING INSTRUCTIONS:

Remove the inner pad from the brace. Rinse the inner pad in water not exceeding 35° with mild detergent for about 5 minutes. Line dry in the shade. Keep out of full sun.

THE FINAL DECISION CONCERNING THE APPLICATION AND CHOICE OF THE EQUIPMENT BELONGS TO A MEDICAL SPECIALIST.

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