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Storage information

The storage location should be dry and at room temperature. The product has to be protected against direct sunlight exposure and humidity.

Cleaning information

Please hang up the bandage dripping wet.

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Textile: Polyester, aluminum.

Disposal

The product can be disposed of with the normal household waste.

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EN

Rhizo Medical thumb orthosis

Thank you very much for selecting a product of the Meyra-Ortopedia Group. The product complies with all legal regulations according to the Medical Products Law.

Please read the directions for use carefully and keep to the advice and recommendations given. Should you have further questions please address the responsible physician, your local specialist dealer or our service hotline.

Product description

The Rhizo medical stabilizing orthosis is a highly effective orthosis on the latest technical level. The material is permeable to air which creates a comfortable microclimate and high wearing comfort. The orthosis stabilizes the thumb metacarpophalangeal joint with a shapeable metal insert at the outer side of the thumb. A steady fit is achieved by means of the two Velcro straps.

Safety instructions

The first application of the Rhizo medical must be carried out by special personnel trained on these kind of products as otherwise the product liability cannot be guaranteed. Basically, when adjusting the Velcro straps it has to be ensured that they are not fastened too tight to avoid congestion.

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Employment according to purpose and indications

The Rhizo medical stabilizing orthosis is used in case of collateral ligament injuries of the thumb metacarpophalangeal joint (skier’s thumb) and irritations in the area of the thumb metacarpophalangeal and thumb saddle joint as well as ligament instabilities.

Application of the bandage

1. Adjust the Velcro fastenings at the thumb to their widest circumference.
2. Slip the orthosis over the thumb while the Velcro strap at the wrist is open. The lateral metal splint runs directly over the thumb joints.
3. Draw the orthosis together at the wrist and fix the Velcro fastenings.
4. Readjust the Velcro fastenings at the phalanx of the thumb fitting the orthosis as accurate as possible. The thumb end phalanx should be movable.
5. Now readjust the wrist band until the orthosis fixes the thumb metacarpophalangeal joint steadily and immobilizes it.

Contraindications

Provided that the product is employed properly contraindications are not know. If the orthosis is too tight pressure symptoms may evolve. In case of the following clinical pictures the product must only applied after agreement with the responsible physician: skin changes or inflammations in the area concerned, scarring, hypersensitivity because of heat production, circulatory disturbance in the hands and arms, lymphatic drainage disturbance.

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Application duration

The wearing period has to be agreed with the responsible physician.

Further information

- The stabilizing orthosis should only be worn under the care of a physician and for the intended clinical pictures.
- The Rhizo medical stabilizing orthosis is intended for the treatment of one patient.
- Improper application excludes the product liability of the manufacturer.
- Should the symptoms aggravate a physician should be consulted immediately.
- A combination with other products has to be agreed with the physician in advance.

Guarantee

The guarantee lies within the framework of the legal regulations.

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