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Cleaning information

Please hang up the bandage dripping wet.

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EN

Achillo Medical Achilles tendon bandage with gel pads

Thank you very much for selecting a product of the Meyra-Ortopedia Group. The product complies with all legal regulations according to the Medical Products Law.

Please read the directions for use carefully and keep to the advice and recommendations given. Should you have further questions please address the responsible physician, your local specialist dealer or our service hotline.

Product description

The Achillo medical Achilles tendon bandage is a highly effective medical bandage on the latest technical level. The anatomically bi-elastic knitted form combined with the gel pads offers a very good fit as well as increased compression, high support and safety in the movement. The Achilles tendon is effectively relieved, the pain eases. The silicone pads improve the blood circulation and support the remission of haematomas and swellings. The enclosed gel heel inlay should always be worn on both sides.

Safety instructions

The first application of the Achillo medical active bandage must be carried out by special personnel trained on these kind of products as otherwise the product liability can not be guaranteed.

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Employment according to purpose and indications

The Achillo medical active bandage is used in the post-treatment of disinsertions of the Achilles tendon as well as Achilles tendon irritations and pain. Any other employment, i.e. departing from the direction, results in an expiration of the product liability.

Application of the bandage

1. Turn the upper edge of the bandage including the silicone cushions down.
2. Now pull the bandage over foot and heel.
3. Then turn the bandage up again and check if the gel pads lay on both sides of the Achilles tendon.
4. Insert the gel heel inlays into the respective shoes.

Contraindications

Provided that the product is employed properly contraindications are not know. If the bandage is too tight pressure symptoms may evolve. In case of the following clinical pictures the product must only applied after agreement with the responsible physician: skin changes or inflammations in the area concerned, scarring, hypersensitivity because of heat production, varicose veins, circulatory disturbance in the legs, lymphatic drainage disturbance.

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Application duration

The bandage should mainly be worn during physical activity as the function only works in this state.

Further information

- The active bandage should only be worn under the care of a physician and for the intended clinical pictures.
- The Achillo medical active bandage is intended for the treatment of one patient.
- Improper application excludes the product liability of the manufacturer.
- Should the symptoms aggravate a physician should be consulted immediately.
- A combination with other products has to be agreed with the physician in advance.

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