

WRIST ORTHOSIS WITH A THUMB HOLD

Catalogue number: MY2761

PRODUCT CHARACTERISTICS:

- stable orthosis for wrist positioning,
- adjustable by means of three hook and loop straps,
- with stable, shapeable aluminium insert,
- fabric permeable to air,
- left and right version,
- colour: black.

USE:

Injuries of the thumb saddle joint, irritations in the wrist area, arthrosis and distorsions, rheumatoid arthritis, instabilities.

COMPOSITION:

PA, POLYESTER, PU FOAM, PP, PE, ALUMINUM

Size	1	2	3	4
Wrist circumference (cm)	13 - 16	16 - 19	19 - 22	22 - 25

